# CareMat

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Bircher Reglomat AG, Wiesengasse 20, CH-8222 Beringen

RoHS 2011/65/EU, GPS2001/95/EG  
CareMat Ax1C, CareMat Bx1C, CareMat Cx1C, CareMat Mx1C

RoHS 2011/65/EU, EMC 2014/30/EU  
CareMat Ax1P, CareMat Bx1P, CareMat Cx1P, CareMat Mx1P

## 11 Kontakt | Seller

**Bircher Reglomat AG**, Wiesengasse 20, CH-8222 Beringen, [www.bircher-reglomat.com](http://www.bircher-reglomat.com) | **Swissdoor ApS**, Stenhuggervej 2, 5471 Søndersø, [www.swissdoor.dk](http://www.swissdoor.dk)

# CareMat

## with replaceable cable

The pressure-sensitive bed mat for dementia patients and persons at risk of falling

### Translation of the original instructions

## 1 Safety instructions

These instructions must be completely read through prior to the commissioning and use of the CareMat. They should be kept safely for future use. Warranty claims will be invalid in the event of any damage to the product, to goods or persons that is due to failure to observe these instructions.

### Instructions for nursing staff

- The CareMat operates purely as a support for the nursing staff, but is not a substitute for them.
- The correct functionality of the CareMat is to be checked each time that it is used.

### Notes for the installer


- The electrical installation work and commissioning of the product may only be carried out by a qualified technician
- The product must not be connected to mains voltage (see "Different Types of Electrical Connections").

## 2 Product description


### 2.2 Areas of application and function

The CareMat was primarily designed for use in retirement and nursing homes, hospitals and psychiatric wards as well as for sheltered residential accommodation. The necessary care and supervision of patients with dementia and of those who are at risk of falling represents a significant increase in care expenditure for care practitioners. The CareMat is the ideal support: Risks of falling and the dangers which arise from wandering or loss of orientation can be significantly reduced.


The CareMat, as a pressure-sensitive sensor, is installed directly in front of the carepatient's or sick person's bed and is connected to the internal call system. As soon as the person leaves the bed and steps onto the CareMat this switches on and sends out an alarm signal to the call system.

**CareMat A01C / A01P**  
**CareMat A11C / A11P**  
1100 x 700 mm

**CareMat B01C / B01P**  
**CareMat B11C / B11P**  
700 x 400 mm


**CareMat C01C / C01P**  
**CareMat C11C / C11P**  
1100 x 700 mm,  
semi-circular


- ① Replaceable cable for connecting to the call system
- ② Open cable end
- ③ Housing with terminals or universal circuit-board
- ④ Chamfered, wheelchair-accessible edges

- ⑤ Non-slip surface with stud texture design
- ⑥ Carrying handle
- ⑦ Active switching zone


## 3 Day to day operation

Position the CareMat centrally in front of the bed, that avoiding the mat is impossible. Route the cable under the bed in order to avoid the risk of tripping over it. Insert the plug into the correct plug socket device.

If the CareMat should not be needed, it can be pushed beneath the bed without the cable connection to the call system having to be disconnected.

Remove CareMat before moving the bed. This will avoid any unnecessary damage to the mat or connecting cable.

## 4 Cleaning and disinfecting


- Disconnect the connection to the call system prior to any cleaning.
- CareMat can be cleaned using commercially available disinfecting and cleaning agents with max. 70 vol% alcohol. No corrosive agents should be used.
- CareMat is neither autoclavable nor can it be disinfected using steam.

## 5 Transport and storage

- During transport and storage any significant bending of the CareMat should be avoided.
- Store the CareMat in a dry, clean and flat location.
- Do not place any heavy loads on the CareMat and do not stack more than 20 packed mats on top of each other.

## 6 Different electrical connections – with terminals circuit-board

The CareMat xx1C (C = cable) includes a terminal circuit-board in the housing.


The correct plug for connecting to the alarm system and its terminal assignment are to be obtained from the relative system manufacturer.

## 7 Different electrical connections – with universal circuit-board

The CareMat xx1P (P = circuit-board) includes a universal circuit-board in the housing.

The universal circuit-board can be adapted to your needs by altering the switch positions. The wiring on the universal circuit-board must not be removed. The connections are protected against wrong polarity of input voltage.


### Variable resistance R

The universal circuit-board has two connector sockets for commercially-available wired resistors. A call resistor that is specifically for the call system for monitoring can be installed in this. The connecting wiring should have a diameter of between 0.38 and 0.56 mm. The socket spacing should be 14 mm.

### NC contact

The NC (= normally closed) contact option that can be selected on the universal PCB is a transistor circuit that operates according to the quiescent current principle. If the CareMat is not actuated, a quiescent current flows through the NC contact as a result of the externally applied voltage. If the resident /

patient actuates the CareMat by stepping on it, this NC contact is opened. This can be detected by a nurse call system operating according to the quiescent current principle that monitors the current flow.

⚠ After the electrical installation, the installer must perform a functional test.

## 8 Maintenance and fault finding

The CareMat is maintenance-free. If the cable should however be pulled out. It can simply be replaced: open housing, replace and connect up cable. The cable ties around the connecting cable must be stripped using a cable tie tool to avoid any slipping. Close the cover and lightly tighten the screws. If despite this, the function is not being provided, the following checks should be made:

- Is the call system operating and is the CareMat correctly connected?
- Is the plug, the plug connection or the cable damaged?

If the points mentioned above have been checked out, but the functions of the CareMat are, despite this, still not being provided, you should consult a qualified expert.

## 9 Disposal

The CareMat can be disposed of in the household refuse.

## 10 Technical specifications

<b>Article</b>	CareMat Ax1C, Ax1P / Bx1C, Bx1P / Cx1C, Cx1P	<b>Cable</b>	Length 3 m, diameter 5.2 ± 0.2 mm
<b>Dimensions (LxWxH)</b>	1100x700x9 mm / 700x400x9 mm / 1100x700x9 mm, height with housing 15 mm	<b>Type of cable</b>	4 pole, LiYYv 4 x 0.34 mm <sup>2</sup>
<b>Weight [kg]</b>	8.2 kg / 3.2 kg / 6.5 kg	<b>Signal output</b>	Mechanical NO contact, transistor NC contact, resistor wiring
<b>Material</b>	Polyurethane (PU)	<b>Electrical load capacity</b>	max. 25 mA at 3 – 24 VDC
<b>Colour</b>	Mat: Careyellow, Caregrey; Cable: white	<b>Type of protection</b>	IP54 without cable end / connecting plug
<b>Surface</b>	Stud texture design	<b>Ambient temperature</b>	0 up to +55 °C
<b>Weight of persons</b>	At least 10 kg		
<b>Chemical resistance</b>	Water (room temperature), body fluids, disinfectant, max. 70 vol% alcohol		<i>Subject to technical modifications</i>

## 10 EC-Declaration of Conformity

Manufacturer:	Bircher Reglomat AG, Wiesengasse 20, CH-8222 Beringen
Following directives have been observed:	RoHS 2011/65/EU, GPS2001/95/EG
Product variants:	CareMat Ax1C, CareMat Bx1C, CareMat Cx1C, CareMat Mx1C
Following directives have been observed:	RoHS 2011/65/EU, EMC 2014/30/EU
Product variants:	CareMat Ax1P, CareMat Bx1P, CareMat Cx1P, CareMat Mx1P

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