

## SAFETY AT EVERY STEP OF THE WAY

**CareMat<sup>®</sup>** THE PRESSURE-SENSITIVE MAT FOR SUPERVISION OF PEOPLE WITH DEMENTIA AND DISORIENTATION

### APPLICATION AND USE IN PATIENT CARE

The CareMat<sup>®</sup> is a pressure-sensitive detection mat, which is connected to the in-house nurse call system. As soon as a person steps onto the mat, a signal is triggered via the nurse call system.



Application areas	<ul> <li>Retirement and nursing homes, hospitals, psychiatric wards and assisted-living facilities</li> </ul>
Installation	<ul> <li>Easy plug-in connection to the nurse call system</li> <li>Can be used in conjunction with the call button</li> </ul>
Features	<ul> <li>Reliable alarm sensor</li> <li>Non-slip material</li> <li>Hygenic surface which is easy to clean</li> <li>Carrying handles for easy handling</li> <li>Patented housing for hard-wired or wireless version</li> <li>Robust and durable construction</li> </ul>
•	<ul> <li>Prevention of falls due to quick assistance</li> <li>Reduces nightly rounds</li> <li>Lowers subsequent health-care costs</li> </ul>
-	<ul><li>Total freedom of movement</li><li>Feeling of safety</li></ul>

### DAILY USE IN A LONG-TERM CARE FACILITY

0.0

Caring for residents with dementia means extra effort for the nursing staff. The CareMat<sup>®</sup> provides ideal support: The risk of falls and the dangers, which can arise from patients with disorientation can be considerably reduced.





Residents with dementia are often restless and disoriented and therefore often try to get up without assistance. This is where CareMat<sup>®</sup> is invaluably.

When a patient steps onto the CareMat<sup>®</sup>, it triggers a nurse call that is relayed to a nurse, who can immediately come to the patient's aid.

The nurse assists the resident and is able to prevent a possible fall by having reacted quickly. This also prevents the resident from leaving the premises without supervision.

CareMat<sup>®</sup> – the daily helper for nursing staff and residents.

### BIRCHER Reglomat

TECHNICAL DATA	CareMat <sup>®</sup> is available in a hard-wired or wireless version for placing beside a bed. The cable respectively the transmitter battery are easily replaceable.	
CareMat A	1100 x 700 x 9 mm, 8.2 kg Housing dimensions 56 x 38 x 10 mm	CareMa
Material	Polyurethane	
Surface	Active grip structure	
Colour	Yellow	
	<ul> <li>Hard-wired: NO or NC contact, resistor circuitry Wireless: Manufacturer specific wireless protocol</li> </ul>	
Connection	<ul> <li>Will be converted in accordance with your nurse call system by your sales representative</li> </ul>	CareMa

# 9

CareMat A01C: Hard-wired



CareMat A01T: Wireless

## A RELIABLE PARTNER, TODAY AND TOMORROW

Bircher Reglomat AG is an internationally leading manufacturer of pressure-sensitive and contact-free detector systems headquartered in Beringen, Switzerland.

The company has been certified in acc. with ISO 9001 and ISO 14001.

CareMat<sup>®</sup> is a quality product and is manufactured in Switzerland.



### Your CareMat® manufacturer

Bircher Reglomat AG Wiesengasse 20 CH-8222 Beringen

Phone +41 (0)52 687 11 11 Fax +41 (0)52 687 11 12 www.bircher-reglomat.com info@bircher.com

### Your CareMat® sales representative